

Lascar EL-P-VACX Hochprzise Glykol-Thermistorsonde, unkalibriert



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Hochprzise, unkalibrierte, 3 Meter lange Thermistorsonde in einer versiegelten, mit Glykol gefüllten Flasche zur Impfstoffberwachung. Zur Verwendung mit dem EL-WiFi-VACX Datenlogger

[\(XE-4058-001\)](#)

Technische Daten

Specifications

Product Code	XE-4082-001
Model	EL-P-VACX
Type	High accuracy Glycol thermistor probe, uncalibrated
Accuracy	+/-0.1°C / +/-0.2°F
Cable Length	3 metres
Measuring range	-40 to +60°C / -40 to + 140°F