


## PCB-Montage-Thermoelement-Stecker FMTC-CU-PCB Typ Cu



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PCB-Montage-Thermoelement-Steckdose FMTC-CU-PCB Typ CU**

Speziell fr Leiterplatten entwickelt, massive Pins fr direkte Platinmontage und Bereitstellung fr CJC-Sensoren

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1379-001
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	PCB socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Brass
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector