



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Typ K Freiliegendes Thermoelement mit geschweifter Spitze 1/0,3 mm PFA Twin Twisted Cable mit nackten Tails 10 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bare-Tails-Terminierung.

Geschweite, freiliegende Verbindung für schnelle Reaktionszeiten. Langlebige Konstruktion mit PFA-Isolierung für Beständigkeit gegen Le, Suren, Flüssigkeiten und gute mechanische Festigkeit und Flexibilität.

Specifications

Product Code	XE-0527-001
General Description	Welded, exposed junction for fast response times. Durable construction with PFA insulation for resistance to oils, acid, fluids and good mechanical strength and flexibility.
Manufacturing Code	Z2-K-10.0-1/0.3-T-JIS
Sensor Type	Exposed Junction
Thermocouple Type	JIS Type K
Cable Length	10m (32.8ft)
Cable Type	PFA Twin Twisted
Core / Strands	1/0.3mm
Cable Termination	Bare Tails
Max. Temperature	+250°C
Min. Temperature	-75°C
Series	Exposed Junctions