

**IEC Typ K 11,8 kg Pull Power (Hufeisen)
Magnet-Thermoelement, PFA-isoliertes Kabel mit
Edelstahl-bergeflecht mit Abschluss im
Miniaturstecker - 1,5 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0075-001
General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	1.5M
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC