



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI Feindraht-Scheiben-Thermoelement, PFA mit Miniaturstecker - Typ K mit Miniatursteckeranschluss - 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0569-001
General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc, 6.35mm diameter x 0.3mm thick
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	5M
Cable Type	1/0.2mm PFA Teflon® insulated twisted wire
Termination Type	Miniature Plug
Max. Temperature	200°C
Min. Temperature	0°C
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow
Accessories	Supplied with high temperature adhesive kapton strips (13x18mm x 5) for quick and easy attachment to surface/process
Series	Copper Disk Thermocouple