



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI Silikonkautschuk-Patch-Thermoelement, PFA Twin Twisted - Typ K mit Miniatur-Steckeranschluss - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0574-001
General Description	Resistant to many oils & chemicals
Sensor Type	Thermocouple junction encapsulated in a silicone rubber patch
Thermocouple Type	K
Cable Length	2M
Cable Type	PFA insulated twin-twist wire
Core / Strands	1/0.2mm
Termination Type	Miniature Plug
Max. Temperature	150°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow
Dimensions	Length 40mm x Width 13mm x Height 5mm Rubber patch
Series	Rubber Patch Thermocouple