

**ANSI Typ K 0,7 kg
Zugknopf-Magnetthermoelement, PFA isoliert mit
Miniaturstecker - 2M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Technische Daten**Specifications**

Product Code	XE-0583-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	1/2" diameter button type, 0.7kg Pull
Cable Length	2M
Cable Type	PFA insulated twin-twist wires
Cable Termination	Miniature Fitted Plug
Max. Temperature	+250°C (plug 200°C)
Standards Met	ANSI MC96.1
Colour	+ Yellow / - Red
Series	Magnet Thermocouple