



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI Bolt Thermoelement, Glasfaser Edelstahl ber geflochtenem Kabel - Typ K - M8 Gewinde - 13mm Schraube



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0584-001
Thermocouple Type	K
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2M
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	250°C
Standards Met	ANSI MC96.1
Bolt thread/pitch	M8 x 1mm
Colour	Positive: Red / Negative: Yellow
Length	13mm Bolt
Series	Bolt Thermocouple