

**ANSI Typ K 1,9 kg
Zugknopf-Magnet-Thermoelement, PFA-isoliertes
Kabel mit bergeflochtenem Edelstahlgeflecht mit
Anschluss im Standardstecker - 1 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0622-001
General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height
Cable Length	1M
Cable Type	Glassfibre insulated lead with stainless steel overbraid 7/0.2mm
Termination Type	Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI
Series	Magnet Thermocouple