

ANSI Krokodilklemmen-Thermoelement mit berflochtenem Glasfaser-Edelstahlkabel - Typ K - 3M Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0629-001
General Description	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Cable Length	3M
Cable Type	Stainless steel over-braided PFA Teflon®
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1
Series	Crocodile Clip Thermocouple