

**ANSI gefertigtes geerdetes Thermoelement aus  
Edelstahlrohr, berflochtenes  
Glasfaser-Edelstahlkabel - Typ K - 4,76 mm  
Durchmesser - 100 mm Länge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                              |
|------------------------------|----------------------------------------------|
| <b>Product Code</b>          | XE-0642-001                                  |
| <b>General Description</b>   | 316 stainless steel tube in various sizes    |
| <b>Thermocouple Type</b>     | K                                            |
| <b>Thermocouple Junction</b> | Grounded at tip for fast response            |
| <b>Cable Length</b>          | 2M                                           |
| <b>Cable Type</b>            | Glassfibre stainless steel overbraided cable |
| <b>Cable Termination</b>     | Bare tails                                   |
| <b>Max. Temperature</b>      | 350°C                                        |
| <b>Min. Temperature</b>      | -60°C                                        |
| <b>Standards Met</b>         | ANSI MC96.1                                  |
| <b>Series</b>                | Fabricated Tube Thermocouple                 |