



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **ANSI-Rohroberflächen-Thermoelement, Glasfaser-Edelstahl-Geflecht - Typ K für Rohrdurchmesser von 22-30 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Leistungsbeschreibung

**Specifications**

<b>Product Code</b>	XE-0655-001
<b>General Description</b>	Stainless steel pipe clamp sensor with adjustable diameter
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded tip for fast response
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre stainless steel braided extension cable
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare Tails
<b>Max. Temperature</b>	+350°C
<b>Standards Met</b>	ANSI MC96.1
<b>To Fit Pipe Diameters</b>	22-30mm
<b>Series</b>	Pipe Surface Thermocouple