



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI-Rohroberflächen-Thermoelement, Glasfaser-Edelstahl-Geflecht - Typ K für Rohrdurchmesser von 40-55 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0657-001
General Description	Stainless steel pipe clamp sensor with adjustable diameter
Thermocouple Type	K
Thermocouple Junction	Grounded tip for fast response
Cable Length	2M
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	ANSI MC96.1
To Fit Pipe Diameters	40-55mm
Series	Pipe Surface Thermocouple