



**LABFACILITY**  
TEMPERATURE & PROCESS TECHNOLOGY

## **ANSI Leaf Thermoelement, glasfaserisoliertes Kabel mit Edelstahlgeflecht, das in blanken Schwänzen endet - 2 m Länge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                                          |
|------------------------------|----------------------------------------------------------|
| <b>Product Code</b>          | XE-0693-001                                              |
| <b>General Description</b>   | Thermocouple mounted on a stainless steel shim           |
| <b>Thermocouple Type</b>     | K                                                        |
| <b>Thermocouple Junction</b> | Junction grounded at tip for fast response               |
| <b>Cable Length</b>          | 2M                                                       |
| <b>Cable Type</b>            | Glassfibre insulated lead with stainless steel overbraid |
| <b>Core / Strands</b>        | 7/0.2mm                                                  |
| <b>Termination Type</b>      | Bare Tails                                               |
| <b>Max. Temperature</b>      | 350°C                                                    |
| <b>Min. Temperature</b>      | -50°C                                                    |
| <b>Standards Met</b>         | ANSI MC96.1                                              |
| <b>Dimensions</b>            | 25x13mm shim (0.2mm thick)                               |
| <b>Series</b>                | Leaf Thermocouple                                        |