



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Allzweck-Thermoelement-Sonde, glasfaserisoliertes Kabel mit Edelstahl-bergeflecht - Typ K mit Anschluss im Standardstecker - 2M Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0765-001
General Description	Rugged thermocouple with stainless steel sheath
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	ANSI MC96.1
Series	Fabricated Tube Thermocouple