


# Standard-rechteckiger Thermoelement-Steckverbinder Faszienbuchse IS-E-FF Typ E IEC



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard-rechteckiger Thermoelement-Steckverbinder Faszienbuchse IS-E-FF Typ E IEC**

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1637-001
<b>General Description</b>	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
<b>Thermocouple Type</b>	E IEC
<b>Connector Type</b>	Fascia Mounting Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Violet
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS