



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

**ANSI-Unterlegscheiben-Thermoelement,
glasfaserisoliertes Kabel mit Edelstahlgeflecht,
das im Standardstecker endet - 4,0 mm
Unterlegscheibenlochgre**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0929-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/0.2mm
Termination Type	Standard plug termination with cable clamps
Max. Temperature	+350°C
Standards Met	ANSI
Washer	4.0mm clearance hole
Series	Washer Thermocouple