



Standard-rechteckiger Thermoelement-Steckverbinder Faszien-Sockel FSTC-CU-FF Typ Cu



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard rechteckiger Thermoelement-Steckverbinder Faszien-Steckdose FSTC-CU-FF Typ CU
 Nylon-Clip für schnelles, einfaches Befestigen am Panel. Alle Fassadenmontage-Fassungen werden durch ungleich große Pins polarisiert
 Spezifikationen

Specifications

Product Code	XE-1677-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	White
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS