



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC-Dsen-Thermoelement, glasfaserisoliertes Kabel mit Edelstahl-Geflecht - Typ K - Miniatur-Steckeranschluss - 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0882-001
General Description	Thermocouple housed in a stainless steel threaded bolt
Thermocouple Type	K
Thermocouple Junction	Grounded hot junction for fast response, M6x1mm thread
Probe Diameter	4.5mm
Probe Length	150mm
Cable Length	5M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Cable Termination	Miniature Plug
Termination Type	Grounded hot junction for fast response, M6x1mm thread
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC
Connector Type	Plug
Connector Size	Miniature
Dimensions	25x13mm shim (0.2mm thick)
Thread	M6 x 1mm thread