




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **PFA-isoliertes, doppelt verdrehtes IEC-Thermoelement mit freiliegender Sperrschicht und blanken Schwänzen - Typ J - 1/0,2 mm - 1 m Kabellänge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement Typ K, J oder T. Schnell ansprechende, geschweifte, freiliegende Verbindungsstelle, erhltlich in der Kabelspezifikation der Klasse 1 gem BS EN 60584-3:2008 mit blankem Ende.

Ideal fr Test- und Entwicklungsanwendungen

### Leistungsbeschreibung

## Specifications

|                     |                                         |
|---------------------|-----------------------------------------|
| Product Code        | XE-7134-001                             |
| General Description | Fast response, welded exposed junction  |
| Sensor Type         | Exposed Junction                        |
| Thermocouple Type   | J IEC                                   |
| Cable Length        | 10M                                     |
| Cable Type          | Teflon® PFA insulated twin-twisted lead |
| Core / Strands      | 7/0.2mm                                 |
| Cable Termination   | Bare tails                              |
| Max. Temperature    | 250°C                                   |
| Min. Temperature    | -75°C                                   |
| Tolerance           | Class 1                                 |
| Standards Met       | IEC-584-3                               |
| Series              | Exposed Junctions                       |