



**LABFACILITY**  
TEMPERATURE & PROCESS TECHNOLOGY

## IEC-Feindraht-Scheiben-Thermoelement, PTFE mit Miniaturstecker - Typ J - 5M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                                              |
|------------------------------|--------------------------------------------------------------|
| <b>Product Code</b>          | XE-3760-001                                                  |
| <b>General Description</b>   | Exposed junction thermocouple with Copper Disc brazed at tip |
| <b>Sensor Type</b>           | Copper disc                                                  |
| <b>Thermocouple Type</b>     | J                                                            |
| <b>Thermocouple Junction</b> | Grounded measuring junction                                  |
| <b>Cable Length</b>          | 5M                                                           |
| <b>Cable Type</b>            | 1/0.2mm PTFE insulated twisted wire                          |
| <b>Termination Type</b>      | Miniature flat pin plug, colour coded in accordance IEC 584  |
| <b>Standards Met</b>         | IEC                                                          |