



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC-Einschraub-Thermoelement-Temperaturfühler - Typ K - 1/4" NPT - 150 mm Länge x 6 mm Durchmesser



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Einschrauben von Temperatursensoren vom Typ K mit einem Messbereich von -75 bis +250 C und 1 Meter PFA-Flachleitungsdraht.

Für den direkten Einbau in Behälter, Rohre, Tanks oder Kanäle zur Messung der Flüssigkeits- oder Gastemperatur.

Typische Anwendungen

Industrielle Prozesskontrolle

Rohre oder Kanäle

Behälter & Tanks

Leistungsbeschreibung

Specifications

Product Code	XE-7354-001
General Description	Screw in Type K thermocouple
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	150mm
Sheath Diameter (mm)	6mm
Process Connection options	1/4"NPT
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	IEC
Number of Wires	2