



## Thermoelement-Flachstecker Typ K



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

- Crimp-auf-Flachstecker
- Preis pro Packung mit 20 Paar
- Geeignet fr IEC-, ANSI- und JIS-Kalibriertypen.
- Hergestellt aus kompatiblen Thermoelmentmaterial
- Thermoelment-Typ, der auf Klemme gestempelt ist
- Suits Labfacility & andere Sperrklemmen
- Fr Drahtgren 0,5 bis 1,0 mm

Klemmenleiste

## Specifications

<b>Product Code</b>	XE-1982-020
<b>General Description</b>	pack of 20 pairs
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Series</b>	SPADE TERMINALS