


Miniatur-Thermoelement-Stecker AM-J-M Typ J ANSI



Miniature in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Stecker AM-J-M Typ J ANSI

Mit flachen Stiften, die zum Thermoelementtyp passen und ideal für den allgemeinen Gebrauch sind, sind alle Kontakte polarisiert, um eine korrekte Verbindung zu gewährleisten. Diese Stecker können Thermoelementkabel mit einem Durchmesser von bis zu 4 mm aufnehmen.

Technische Daten**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-9060-001 |
| General Description | Ideal for general purpose use. Contacts are polarised to ensure correct connection. |
| Thermocouple Type | J ANSI |
| Connector Type | Plug |
| Connector Size | Miniature |
| Colour | Black |
| Max. Temperature | 220°C |
| Pins | Flat Pins |
| +Positive Contact | Iron |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |
| Series | Miniature Connector |