



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC-gefertigte Rohrthermoelemente mit Crimpdichtung und Anschlussdraht - Typ K - 4,5 mm Manteldurchmesser - 100 mm Mantellänge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typische Anwendungen

Industrielle Prozesskontrolle

Rohre oder Kanäle

fen und Heizung

Universell einsetzbar

Leistungsbeschreibung

Specifications

Product Code	XE-7905-001
General Description	Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	K
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	100mm
Sheath Diameter (mm)	4.5mm
Max. Temperature	250°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC
Number of Wires	2
Series	Fabricated Tube Thermocouple