




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC-gefertigte Rohrthermoelemente mit Crimpdichtung und Anschlussdraht - Typ J - 4,5 mm Manteldurchmesser - 100 mm Mantellänge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typische Anwendungen
Industrielle Prozesskontrolle
Rohre oder Kanle
fen und Heizung
Universell einsetzbar

Leistungsbeschreibung**Specifications**

| | |
|-----------------------------|---|
| Product Code | XE-7917-001 |
| General Description | Fabricated IEC Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead |
| Thermocouple Type | J |
| Cable Length | 1M |
| Cable Type | PFA |
| Termination Type | Stripped Wires |
| Sheath Length (mm) | 100mm |
| Sheath Diameter (mm) | 4.5mm |
| Max. Temperature | 250°C |
| Min. Temperature | -75°C |
| Tolerance | Class 1 |
| Standards Met | IEC |
| Number of Wires | 2 |
| Series | Fabricated Tube Thermocouple |