


Miniatur-Thermoelement-Steckverbinder Inline-Buchse AM-J-FS Typ J ANSI



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Steckverbinder Inline-Buchse AM-J-FS Typ J ANSI

Ideal für den allgemeinen Gebrauch, sind alle Kontakte polarisiert, um eine korrekte Verbindung zu gewährleisten. Diese Stecker können Thermoelementkabel mit einem Durchmesser von bis zu 4 mm aufnehmen.

Technische Daten**Specifications**

Product Code	XE-9062-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	J ANSI
Connector Type	Socket
Connector Size	Miniature
Colour	Black
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS