

ANSI Fabricated Tube Thermoelemente mit Crimpdichtung und Anschlussdraht - Typ K - 6,0 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typische Anwendungen

Industrielle Prozesskontrolle

Rohre oder Kanäle

fen und Heizung

Universell einsetzbar

Leistungsbeschreibung**Specifications**

Product Code	XE-7947-001
General Description	Fabricated Type K Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	K
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	250
Sheath Diameter (mm)	6.0
Max. Temperature	250°C
Min. Temperature	75°C
Tolerance	Class 1
Standards Met	ANSI MC96.1
Number of Wires	2
Series	Fabricated Tube Thermocouple