



## Miniatur Quick Wire Thermokoppel-Stecker AM-E-MQ Typ E ANSI



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Stecker Schnelldrahtstecker AM-E-MQ Typ E ANSI**

Miniaturstecker mit flachen Steckerstiften. Einfache Jab-In-Verbindung fr eine schnelle Kndigung. Einfach den Draht eindrcken und die Schraube festziehen.

**Technische Daten****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9041-001                                   |
| <b>General Description</b> | Easy jab-in connection for rapid termination. |
| <b>Thermocouple Type</b>   | E ANSI  |
| <b>Connector Type</b>      | Quick Wire Plug                               |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Violet  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>Pins</b>                | Flat Pins                                     |
| <b>+Positive Contact</b>   | Nickel Chromium                               |
| <b>-Negative Contact</b>   | Copper Nickel                                 |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |