



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

PTFE-isoliertes doppelt verdrehtes IEC-Offenschaltungs-Thermoelement mit nackten Endenden Typ K 300 mm Kabellänge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement Typ K, J oder T. Schnelle Ansprechanlage, geschweite freiliegende Verbindung, verfügbare nach der Klasse-1-Kabelspezifikation gem BS EN 60584-3:2008 mit Bare-Tails-Abschluss.

Ideal für Test- und Entwicklungsanwendungen.

Technische Daten

Specifications

Product Code	XE-0424-001
General Description	Fast response, welded exposed Junction
Sensor Type	Exposed Junction
Cable Length	300mm
Core / Strands	1/0.2mm
Cable Termination	Bare tails
Max. Temperature	+260°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Series	Exposed Junctions