


# Miniatur Quick Wire Thermoelement-Steckverbinder AM-J-FQ Typ J ANSI



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Quick Wire Connector Thermoelement-Buchse AM-J-FQ Typ J ANSI**

Einfache Einstichverbindung, einfach den Draht eindrücken und die Schraube festziehen. Die Kontakte sind polarisiert, um eine falsche Terminierung zu verhindern.

**Technische Daten****Specifications**

<b>Product Code</b>	XE-9072-001
<b>General Description</b>	Easy jab-in connection, just push in wire & tighten screw.
<b>Thermocouple Type</b>	J ANSI
<b>Connector Type</b>	Quick Wire Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Black
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Iron
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS