

RFX Drahtlose Fleischsonde



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Verfolgen Sie die Kernfleischtemperaturen von Ihrem Gert aus

Verwenden Sie bis zu 50 Sonden mit einem GATEWAY

Zuverlässige drahtlose Verbindung

50 Stunden Akkulaufzeit mit einer 10-minütigen Aufladung

Vollständig wasserdicht nach IP69K

Kostenloses rückführbares Kalibrierzertifikat

Muss mit dem

oder

gekoppelt werden

[RFX Gateway](#)
[RFX Starter Kit](#)

Leistungsbeschreibung

Specifications

Product Code	XE-4201-001
Model	RFX Wireless Meat Probe
Type	Thermistor, 4 Sensors
Accuracy	+/-0.5°C (-10 to 100°C)
Battery	1.5 V AAA with charging block - 175+ full charges per AAA battery. 0-80% in 10 minutes; 0-100% in 20 minutes
warranty	6 Months
Output	Wireless communication: 433Mhz radio frequency
Temperature Range	Internal: -18 to 100°C / External: -18 to 538°C
Standards Met	IP69K
Measuring range	Battery powered RFX GATEWAY: 650m (2,132 ft.), AC Powered RFX GATEWAY: 385m (1,263 ft.) (RFX MEAT to RFX GATEWAY)