



**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

**ANSI  
Einschraub-Thermoelement-Temperaturfühler -  
Typ K - 1/4" BSPP - 150 mm Länge x 3 mm  
Durchmesser**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Einschrauben von Temperatursensoren vom Typ K mit einem Messbereich von -75 bis +250 C und 1 Meter PFA-Flachleitungsdraht.

Für den direkten Einbau in Behälter, Rohre, Tanks oder Kanäle zur Messung der Flüssigkeits- oder Gastemperatur.

Typische Anwendungen

Industrielle Prozesskontrolle

Rohre oder Kanäle

Behälter & Tanks

**Leistungsbeschreibung**

## Specifications

<b>Product Code</b>	XE-7366-001
<b>General Description</b>	Screw in Type K thermocouple
<b>Thermocouple Type</b>	K ANSI
<b>Cable Length</b>	1M
<b>Cable Type</b>	PFA
<b>Termination Type</b>	Stripped Wires
<b>Sheath Length (mm)</b>	150mm
<b>Sheath Diameter (mm)</b>	3mm
<b>Process Connection options</b>	1/4"BSPP
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI
<b>Number of Wires</b>	2