


Miniaturrechteckige Faszien-Steckdose Thermoelement-Steckverbinder AM-N-FF Typ N ANSI



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniaturrechteckiger Thermoelement-Steckverbinder Fassadenbuchse AM-N-FF Typ N ANSI

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

Technische Daten**Specifications**

Product Code	XE-1277-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	N ANSI
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Orange
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS