



Typ K PVC isoliertes flaches Thermoelement-Kabel / Draht (IEC) 13/0,2 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement-Kabel / Draht (IEC) Typ K PVC-isoliertes Flachpaar

Das Kompensierungskabel hñelt der thermoelektrischen Eigenschaft des jeweiligen Thermoelements bei einem begrenzten Umgebungstemperaturbereich, typischerweise 0 bis 80C. Da diese Kabel aus kostengünstigen Materialien gefertigt sind, können bei Anlageninstallationen im Vergleich zum Verlegen eines echten Thermokoppel-Verlängerungskabels Kosteneinsparungen erzielt werden.

Technische Daten
Specifications

General Description	Good general purpose insulation for 'light' environments. Waterproof and very flexible.
Cable Length	10, 25, 50, 100 or 200m
Cable Type	Flat pair construction
Core / Strands	13/0.2mm
Conductor Cross Section 13/0.2 (mm²)	0.408
Outer Colour	Green
+ Colour	Green
- Colour	White
Drain Wire	No
Insulation Material	PVC Insulated
Loop Resistance 13/0.2 (Ohms per combined metre)	4.5
Max. Temperature	105°C
Min. Temperature	-10°C



+ Conductor	Nickel Chromium
- Conductor	Nickel Aluminium
Screen	No
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	K
Tolerance	Class 2