


PCB-Montage-Thermoelement-Stecker Steckdose AM-T-PCB Typ T ANSI Zinnbeschichtete Pins



Miniature PCB Socket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PCB-Montage-Thermoelement-Stecker Buchse AM-T-PCB-TIN Typ T ANSI

Speziell fr Leiterplatten entwickelt, massive Zinnstifte fr direkte Platinenmontage und Bereitstellung fr den CJC-Sensor

Technische Daten

Specifications

Product Code	XE-1230-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	T ANSI
Connector Size	Miniature
Connector Type	PCB socket
Colour	Blue
Max. Temperature	220°C
Pins	Tin Plated
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS