



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Mineral Insulated Thermoelement mit Miniaturstopfen Typ K 1,5 mm x 500 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Technische Daten**Specifications**

Product Code	XE-2908-001-1
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Mineral Insulated
Thermocouple Type	K
Element/hot junction	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
Probe Diameter	1.5mm
Probe Length	500mm
Sheath Material	Stainless steel 310
Termination Type	Miniature flat pin plug, colour coded in accordance JIS
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	JIS
Connector Type	Plug
Connector Size	Miniature
Colour	Blue
Series	Mineral Insulated