



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## JIS Mineralisiertes Thermoelement mit Standard-Stecker Typ K 1,5 mm x 150 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

### Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS) Kalibrierungsstandards rckverfolgbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-Rckverfolgbarkeit ber unsere hauseigene Kalibrierungseinrichtung.

**Bitte beachten Sie, dass Sie, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen, bezglich Lieferinformationen kontaktiert werden.**

### Technische Daten

## Specifications

<b>Product Code</b>	XE-2912-001
<b>General Description</b>	Highly flexible, sheath can be bent/formed to suit many applications and processes
<b>Sensor Type</b>	Mineral Insulated
<b>Thermocouple Type</b>	K
<b>Element/hot junction</b>	Single element, junction insulated from sheath (offers protection against spurious electrical signals)
<b>Probe Diameter</b>	1.5mm
<b>Probe Length</b>	150mm
<b>Sheath Material</b>	Stainless steel 310
<b>Termination Type</b>	Standard round pin plug, colour coded in accordance JIS
<b>Max. Temperature</b>	1100°C
<b>Min. Temperature</b>	-40°C
<b>Standards Met</b>	JIS
<b>Connector Type</b>	Plug



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Connector Size	Standard
Colour	Blue