



Miniatur-Panelmontage-Thermoelement-Anschlus s mit Edelstahlhalterung AM-R/S-SSPF Typ R/S ANSI



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Panelmontage-Thermoelement-Anschluss mit Edelstahlhalterung AM-R/S-SSPF Typ R/S ANSI

Edelstahlhalterung, die einfach am Verbinder befestigt ist. Vollständig geschlossene Anschlüsse mit Kabelanschlussdichtung

Technische Daten**Specifications**

Product Code	XE-1318-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	R/S ANSI
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS