




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Typ K PTFE Flachpaar-Thermoelement-Kabel / Draht IEC - 1/0,315 mm Strang



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typ K PTFE isoliertes flachpaariges Thermoelement-Kabel / Draht (IEC)

PTFE-isoliertes flaches Paar. PTFE besitzt eine gute mechanische Festigkeit und Flexibilität und ist resistent gegen Le, Suren und andere unangenehme Flüssigkeiten.

Technische Daten

Specifications

General Description	Resistant to oils, acids adverse agents and fluids. Good mechanical strength and flexibility.
Cable Length	10m, 25m, 50m, 100m, 200m
Cable Type	Flat pair cable
Core / Strands	1/0.315mm
Conductor Cross Section 1/0.315 (mm²)	0.078
Outer Colour	Green
+ Colour	Green
- Colour	White
Drain Wire	No
Insulation Material	PTFE insulated
Loop Resistance 1/0.315 (Ohms per combined metre)	12.8
Max. Temperature	260°C
Min. Temperature	-75°C
+ Conductor	Nickel Chromium
- Conductor	Nickel Aluminium



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Screen	No
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	K
Tolerance	Class 1