


## Typ K PFA isoliertes flachpaariges Thermoelement-Kabel / Draht (JIS) 1/0,315 mm Strang



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PFA-Teflon-isoliertes Flachpaar. PFA-Teflon besitzt eine gute mechanische Festigkeit und Flexibilität und ist widerstandsfähig gegen Öl, Säuren und andere ungünstige Flüssigkeiten, Gase, Dampf und Druck.

### Technische Daten

## Specifications

<b>General Description</b>	Resistant to oils, acids other adverse agents and fluids. Good mechanical strength and flexibility.
<b>Cable Length</b>	10m, 25m, 50m, 100m, 200m
<b>Cable Type</b>	Single pair of PFA insulated conductors, laid flat with PFA jacket
<b>Core / Strands</b>	1/0.315mm
<b>Outer Colour</b>	Blue
<b>+ Colour</b>	Red
<b>- Colour</b>	White
<b>Insulation Material</b>	PFA Insulated
<b>Loop Resistance 1/0.315 (Ohms per combined metre)</b>	12.8
<b>Max. Temperature</b>	260°C
<b>Min. Temperature</b>	-75°C
<b>+ Conductor</b>	Nickel Chromium
<b>- Conductor</b>	Nickel Aluminium
<b>Screen</b>	No
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Thermocouple Type</b>	K