



PCB-Montage-Thermoelement-Steckdose AM-K-PCB Typ K ANSI



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PCB-Montage-Thermoelement-Steckdose AM-K-PCB Typ K ANSI

Speziell fr Leiterplatten entwickelt, massive Pins fr direkte Platinmontage und Bereitstellung fr CJC-Sensoren

Technische Daten**Specifications**

Product Code	XE-1139-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	K ANSI
Connector Type	PCB socket
Connector Size	Miniature
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector