


# PCB-Montage-Thermoelement-Steckbuchse AM-E-PCB Typ E ANSI



**Miniature PCB Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PCB-Montage-Thermoelement-Steckbuchse AM-E-PCB Typ E ANSI**

Speziell fr Leiterplatten entwickelt, massive Pins fr direkte Platinmontage und Bereitstellung fr CJC-Sensoren

**Technische Daten****Specifications**

<b>Product Code</b>	XE-9039-001
<b>General Description</b>	Versatile, high quality circuit board mounting miniature socket
<b>Thermocouple Type</b>	E ANSI
<b>Connector Type</b>	PCB socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Violet
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS