


Standard-Rechteckige Faszien-Thermoelement-Steckersteckdose AS-N-FF Typ N ANSI



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standardmiger rechthoekiger Thermoelement-Steckverbinder Faszien-Sockel AS-N-FF Typ N ANSI

Nylonclip fr schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich groe Pins polarisiert

Technische Daten**Specifications**

Product Code	XE-1577-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	N ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Orange
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS