


Standard-Thermoelement-Steckverbinder mit Edelsthalterung AS-N-SSPF Typ N ANSI



Standard in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Steckverbinder mit Edelstahlhalterung AS-N-SSPF Typ N ANSI

Fr den Einsatz in industrielleren Prozessen gedacht, kann aber in allen Anwendungen nach Wunsch eingesetzt werden

Technische Daten**Specifications**

Product Code	XE-1578-001
General Description	Intended for use in more industrial processes but can be used in all applications as desired
Thermocouple Type	N ANSI
Connector Type	Socket with SS Bracket
Connector Size	Standard
Colour	Orange
Max. Temperature	220°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS