



Miniatur-Thermoelement-Anschluss In-Line-Buchse JM-K-FS Typ K JIS



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Anschluss In-Line-Buchse JM-K-FS Typ K JIS

Ideal for general use, all contacts are polarised, to ensure correct connection.
These connectors can accept thermocouple cables with a diameter of up to 4 mm.

Technische Daten**Specifications**

Product Code	XE-9152-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	K JIS
Connector Type	Socket
Connector Size	Miniature
Colour	Blue
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS