


Miniatur-Thermoelement-Steckverbinder In-Line-Buchse JM-T-FS Typ T JIS



Miniature in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Steckverbinder In-Line-Buchse JM-T-FS Typ T JIS

Ideal fr den allgemeinen Gebrauch, sind alle Kontakte polarisiert, um eine korrekte Verbindung zu gewhrleisten.
Diese Stecker knnen Thermoelementkabel mit einem Durchmesser von bis zu 4 mm aufnehmen.

Technische Daten**Specifications**

Product Code	XE-9252-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
Thermocouple Type	T JIS
Connector Type	Socket
Connector Size	Miniature
Colour	Brown
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector