



Miniatur Quick Wire Thermoelement Steckdose JM-B-FQ Typ B JIS



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur Quick Wire Connector Thermoelement-Buchse JM-B-FQ Typ B JIS

Einfache Einstichverbindung, einfach den Draht eindrücken und die Schraube festziehen. Die Kontakte sind polarisiert, um eine falsche Terminierung zu verhindern.

Technische Daten**Specifications**

Product Code	XE-9223-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	B JIS
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Grey
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS