


Miniaturrechteckige Faszien-Steckdose Thermoelement-Anschluss JM-E-FF Typ E JIS



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniaturrechteckiger Thermoelement-Steckverbinder Faszien-Sockel JM-E-FF Typ E JIS

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

Technische Daten**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9237-001 |
| General Description | All fascia mounting sockets are polarised by means of unequally sized pins |
| Thermocouple Type | E JIS |
| Connector Type | Fascia mounting socket |
| Connector Size | Miniature |
| Colour | Violet |
| Max. Temperature | 120°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |