


# Miniatur-Panelmontage-Thermoelement-Anschlus s mit Edelstahlhalterung JM-K-SSPF Typ K JIS



**Miniature in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Panelmontage-Thermoelement-Anschluss mit Edelstahlhalterung JM-K-SSPF Typ K JIS**

Edelstahlhalterung, die einfach am Verbinder befestigt ist. Vollständig geschlossene Anschlüsse mit Kabelanschlussdichtung

**Technische Daten****Specifications**

<b>Product Code</b>	XE-9158-001
<b>General Description</b>	Fully enclosed terminals with cable port seal
<b>Thermocouple Type</b>	K JIS
<b>Connector Type</b>	Socket with SS Bracket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Blue
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS